

Crocus <http://www.crocusscientific.com>
Scientific Trading Tel: +92 42 35786141
173-B Raja Center Main Market Gulberg II Lahore PAKISTAN
E-mail: crocus.scientific@gmail.com

Activity Wheel CST-AW02

Crocus's The Rodent Activity Wheel represents a very simple and clever way to register animal voluntary physical activity in its home cage environment. The use of this high throughput tool is particularly relevant for research involving circadian rhythms, phenotyping and drug testing. The animals are housed individually in the home cages equipped with the running wheel. Two models of rodent activity wheel are available with external or internal wheel. The external wheel option is provided to preserve animal living space.

KEY FEATURES

Easy way to quantify rodent voluntary exercise in their home cage environment
External or internal wheel options
For rat, mice and hamsters
Ideal for high throughput experiments
**Manufacturer: Crocus Scientific Trading
Lahore Pakistan**



The Mouse Tent CST-MT02 is made from tinted polycarbonate and provides an ideal shelter for mothers to build their nests. Its unique design increases the floor area and offers life saving refuge during cage flooding events by providing an opportunity for escape to a higher level where mice can stay warm and dry. • Increases cage floor area • Autoclavable red polycarbonate • Cagewasher safe • Excellent quality yet highly economical • Easy to clean, store & stack.

Specification:

Floor area 100 x 110 mm Height 75mm

