

# Crocus Scientific Trading

## HOUSING

**crocus scientific trading** Manufacturer and Supplier Of Mice & Rats Cages for animal welfare and meeting research requirements in Universities and Research Centers in Pakistan. **crocus** is the ideal partner at your side for your housing needs. we provide solutions for any conventional housing requirement designed for husbandry to meet hygiene needs, occupational health and safety requirements, everything undergoes the top quality check to deliver the reliability you expect from us for a long-lasting product.


## CONVENTIONAL CAGES

“**crocus**” cage bodies are designed to facilitate animal welfare, meet research requirements and minimize experimental variables. the following cage body selection is grouped by species. transparent plastic with exceptional impact strength and heat resistance.it can be autoclaved at 121°C/250°F

### Body:-

- *Autoclave able.*
- *Fully Transparent Body.*
- *Food Grade Quality Material polycarbonate Material.*
- *Excellent stake ability for ease of handling and storage.*
- *Chrome plated/S.S Top grill with duly spot welded and finished joints.*
- *Provision of fixing water bottle.*
- *Provision for feed pallets.*

Manufactured By : **Crocus Scientific Trading**



### Accessories:-

## BOTTLE SELECTION

“**crocus**” bottles have a long lasting silicone ring enabling easy placing and removal of caps; the ring doesn't have to be removed for cleaning when bottles are autoclaved and features three protruding lips which prevent water from leaking. Bottles are made of clear-glass Polycarbonate (PC), FDA requirements and Autoclavable at 121°C /250°F for 20 minutes. Water Bottles Made Of Glass . Glass Bottles Are Graduated From 250ml To 500ml, And Slipper Tube Is Provided In Stainless Steel. leak-proof bottles with silicon rings feature polished with cone shaped, design that maximize water availability for animals;



# Crocus Scientific Trading

## Stainless Steel Feeding Needles

**“crocus”** stainless steel feeding needles meet the needs of researchers that prefer traditional reusable metal tubes. The needles are comprised of a stainless steel tube and ball and a nickel-plated brass hub. The sizes range to cover rodents ranging in size from mouse pups to large rats. The needles are sold straight, but are easily bent for optimized use.



## Trolley for Mice & Rats Cages

**crocus scientific trading** discover the widest range of conventional housing solutions promoting animal welfare and meeting research requirements. CST-T1 Racks/Trolleys for holding animal Cages either mice or rats are available in hanging types. In hanging type, cage is kept in the rack.

### Specifications:-

Size of trolley:

Height=63.5”approx

Width=69”approx

Depth=11”approx.

Cages Capacity= 25 cages

Number of cages in /Row =5

Total shelves= 5



\*Crocus Scientific Trading reserves the right to change in design without notice.

S.No	Item Name	Size	Models No
1	Rats Cage Breeding ss304	419x280x178	CST-RCB-304-1
2	PREMIUM BREEDING SS304 6mm	419x280x178	CST-PRB-6-1
3	Rats Cage Breeding chrome	419x280x178	CST-RCB-CH-1
4	Rats Cage Normal Cage chrome	419x280x178	CST-RCN-1
5	Rats Cage normal SS304	419x280x178	CST-RCN-304-1
6	Mice Cage Breeding ss304	290x220x140	CST-MCB-304-1
7	PREMIUM SS304 6MM	290x220x140	CST-PMB-1
8	Mice Cage Breeding chrome	290x220x140	CST-MCB-CH-1
9	Mice Cage Normal Cage chrome	290x220x140	CST-MCN-CH-1
10	Mice Cage normal SS304	290x220x140	CST-MCN-304
11	Single Mice Cage SS304	230x180x125	CST-M-304-1
12	Single Mice Cage Chrome	230x180x125	CST-M—1
13*	Extra Large mice and Rats Cages	510x330x178	CST-XMR-304

# Crocus Scientific Trading

## NEW ARRIVAL

### Extra Large Mice and Rats Cages

