

Crocus <http://www.crocusscientific.com>

Scientific Trading

Tel: +92 42 35786141

173-B Raja Center Main Market Gulberg II Lahore PAKISTAN

E-mail: crocus.scientific@gmail.com

Automated Y Maze

Description

Crocus Y maze is an enclosed apparatus in the form of a Y placed horizontally, similar to the Crocus T maze. Animals usually start from a stem of the Y Maze (arbitrarily chosen) and allowed to spontaneously alternate between the arms. The test relies on either spontaneous alternation or rewarded alternation (1). The Crocus Maze automated Y maze uses automated detection so that the mouse or rat can be held in the center area with a freeze time, allowing for many trials without human intervention, maximizing ease of use and automated experiments.

Specification

Mouse

- 30cm arm length
- 5cm arm width
- 15cm arm height
- Material: Wood & Acrylic

