

# Crocus

## Mice Rat Pole Apparatus

### Mice Rat Pole Apparatus

The Rat Pole Climbing Apparatus (also known as Cook's Pole Climbing Apparatus) is a cornerstone in laboratory research for evaluating the impact of psychotropic and antipsychotic drugs on rodents. This device simulates a conditioned avoidance response, helping researchers measure learning, memory, and motor coordination in rats and mice. Special grade back light climbing pole for stability and grip.

Model . CST03

Dimensions 16"x16"x16"

Material Acrylic

**Crocus scientific trading**

+92 3228074166

[crocus.scientific@gmail.com](mailto:crocus.scientific@gmail.com)

[www.crocusscientific.com](http://www.crocusscientific.com)

