


**Equipped with its
Unique Feature
Set**



Hot Plate Analgesia Meter

Specifications:

- Surface temperature range: 25C to 80C precisely controlled to within 0.2C
- Built-in electronic timer
- Digital set point
- Optional Data Transfer software
- Built-in digital thermometer
- Start/stop foot switch control
- Temperature range: 25C to 80C
- Temperature Accuracy: 0.1C
- Heated Surface: 3/4" Anodized Aluminum provides a high temperature stability and even surface distribution
- Accuracy of Timer: 0.01 seconds
- Size of heated Surface: 10"x10" x 0.75" (254mm x 245mm x 19mm) - Rs-232 , Wifi transfer of data
- Thermal printer (optional)

Features:

Plate temperature can be held at a set point between 5°C above ambient to 80°C ($\pm 0.1^\circ\text{C}$)
Instrument has provision to add max sequence numbers, auto experiment number and animal sex.
Precise and reproducible results
Cut off timer to avoid damage to animal
PC Connectivity with software
Software with calibration and report generation feature
Software has provision to add title, comment, graphical presentation and data can be converted in excel or pdf

Crocus Scientific Trading

WhatsApp: 0333-4224578, info@cst.com