

Crocus

## Convulsimeter

### Electro Convulsimeter

Convulsimeter (Electro Convulsometer)

Output: Supramaximal electroshock (MES) for rodents.

Output Current: Adjustable 0-150 mA.

Output Duration: 0.1-1.0 sec (adjustable).

Electrodes: Corneal and ear-clip type, with flexible leads.

Digital display for current and time.

Power: 220-240 V.

**Model . CST-DEC-03**

**Crocus scientific trading**

+92 3228074166

[crocus.scientific@gmail.com](mailto:crocus.scientific@gmail.com)

[www.crocusscientific.com](http://www.crocusscientific.com)

